

TubeSeq Ultimate

55 Sample Submission Guide for the TubeSeq Ultimate Service

DNA & PRIMER CONCENTRATION & VOLUME

Plasmid DNA	<30 kbp	50- 100 ng/µl	15 µl
BAC/PAC/ Cosmid DNA	30- 200 kbp	200- 1000 ng/µl	15 µl
Purified PCR Products	150-300 bp 300-1000 bp 1000-3000 bp	1 ng/µl 5 ng/µl 10 ng/µl	15 µl
Unpurified PCR Products	150-300 bp 300-1000 bp 1000-3000 bp	4 ng/µl 10 ng/µl 20 ng/µl	15 µl
Enclosed Primer	16-25 bases	10 pmol/µl	15 µl*


*Sufficient for up to 200 reactions.; $30 \ \mu$ l are needed for up to 350 reactions. For more reactions, the primer must be entered several times and with a new name (e.g. indexed)

OPTIMUM PRIMER CONDITIONS

- The optimum primer length is 16-25 bases
- Primer melting temperature (Tm) should be 50-62° C
- The GC content of the primer should be 35-60 %
- Ideally one G or C should be located at the 3' primer end
- The number of 3' Gs or Cs should not exceed 2 Gs or Cs
- If possible, avoid >3 identical bases in a row in the sequence
- Primers must not contain phosphorylation or fluorescent dyes
- Primer must be disolved in double distilled water or 5 mM Tris-HCl

Please send your samples to:

Eurofins Genomics Sequencing Lab Anzingerstraße 7a / DE - 85560 Ebersberg



THE DNA UNIVERSE

Sample Preparation & Submission.







3. SUBMIT SAMPLES

options to send



4. DOWNLOAD RESULTS You will receive an email once your data is ready to your samples via mail download.

> Data will be available in your online account.

Please ensure to send your samples directly to our sequencing lab in Ebersberg by using one of our shipping options!

1. DEFINE REACTIONS

Login to your Eurofins Genomics account. Choose "TubeSeg Ultimate" in the order menu.

Define your reactions on the order page and allocate the barcodes to your reactions.

2. LABEL TUBES Affix the barcode on your sample tube. Use our Primer Barcodes for enclosed primers.

to seal your tubes!

Use one of our shipping

or courier. Make sure to pack your Please don't use parafilm samples safely.

Additional Service & Information

>>> If you need primer barcodes for your enclosed primers, get our Primer Barcodes here

>> Shipping options for sending samples to our sequencing lab in Ebersberg can be found here

>>> Do you need boxes for safe shipment of your sample tubes? Visit us here



USE OUR FREE SEQUENCING PRIMER DESIGN TOOL TO DESIGN THE OPTIMUM **SEQUENCING PRIMERS**

If you search for the perfect sequencing primer, you will love our SeqPrimer and Standard Primer. Test them with a 25,- Euro voucher (or equivalent national currency). Simply use the promotion code "OLIGOPLANET" during check-out.

