

Sequencing order sheet



Contact details

Customer:	<input type="text"/>	For internal use:	<input type="text"/>
PO number (optional):	<input type="text"/>		<input type="text"/>
Company/Institute:	<input type="text"/>		<input type="text"/>
Department:	<input type="text"/>		<input type="text"/>
Address/Street:	<input type="text"/>		
Post code/City:	<input type="text"/>		
Phone:	<input type="text"/>	Fax:	<input type="text"/>
E-Mail:	<input type="text"/>	Date:	<input type="text"/>

Primer

concentration 10 μ M, in H₂O or TE

No.	Name and sequence:	T _M

Standard primer are available. Please do not include them into the order.

Template

Plasmid: at least 800ng DNA per run, in H₂O. concentration at least 150-200 ng/ μ l, inserts to max. 5 kbp

PCR-product: 100ng DNA per 500 bp length per run

plasmid dissolved in
 PCR-product dissolved in

App. length:



No.	template (max 15 characters)	conc. template [μ g / μ l]	volume template [μ l]	primer	Add. information (i.e. high GC, repeats, vector name, etc.)

Order

Data delivery via: Email SSL-encrypted web-download

To: Nava Harish (Navah@synteza.com) Synteza Bioscience Ltd., Hartom 10, Har Hotzvim Jerusalem
 Tel. 02-5867138 ext. 130 or 02-5867138 ext. 105
 Technical support: Navah@synteza.com or Noah@synteza.com